

2X NeoTaq Master Mix

Features

- ❖ Robust amplification with high yields of PCR product.
- ❖ Minimal optimization due to uniquely formulated buffer.
- ❖ High functional activity.
- ❖ Enhanced efficiency, specificity, and sensitivity.
- ❖ Amplification of long targets up to 6.4kb from genomic DNA.

Applications

- ❖ End point PCR with low copy targets.
- ❖ Amplification from multiple template sources.
- ❖ Multiplex primer reaction.
- ❖ High throughput PCR Procedures.

Quality Control Assays

- ❖ Nuclease assays: No detectable endonuclease, exonuclease and RNase activity.
- ❖ *E. coli* Host contamination: No *E. coli* DNA contamination was detected in qPCR with specific primers targeting 16S rRNA gene.
- Functional Assay: NeoTaq master tested extensively for its reproducible performance in critical PCR amplifications.


Product Description

2X NeoTaq master mix is a ready-to-use reaction mix optimized for routine PCR applications. It contains high quality NeoTaq DNA polymerase, dNTPs, MgCl₂, enhancers and stabilizers. All the reagents required for PCR (except template and primer) is at optimal concentrations without the need of additional optimization step. 2X NeoTaq master mix efficiently amplifies a wide range of DNA templates for most of the PCR applications.

Characterization Studies

➤ Consistent yields of PCR product

E. coli genomic DNA were amplified using 2x NeoTaq master mix and equal volumes of PCR products were analyzed on agarose gel. Lane 1: Ladder; Lane 2-6: Respective PCR product.



➤ Tolerance of variable magnesium ions.

PCR amplifications was performed using 2X NeoTaq master mix with different concentrations of magnesium ions using *E. coli* genomic DNA. PCR products were analyzed on agarose gel.



Lane 1: Ladder; Lane 2: 2mM MgCl₂; Lane 3: 3mM MgCl₂ Lane 4: 4mM MgCl₂

Advantages

- ❖ Reduce the risk of contamination from multiple pipetting steps.
- ❖ Decreases assay set up time.
- ❖ Consistent reaction-to-reaction performance.

ISO9001:2015




ISO13485:2016



➤ Specific amplification with different sources


Three differently sized products from genomic DNA were amplified from multiple template sources using inhouse and commercial 2x PCR master mix. Only 10% of each reaction was loaded on the gel. M: Marker.



Lane 1: Ladder; Lane 2: Human gene(230bp); Lane 3: Plant gene(300bp); Lane 4: E. coli gene(1kb)

Benchmarking studies: *Thermal inactivation (97.5°C)*

2X NeoTaq Master mix is tested for its thermal inactivation property and found to be best among the 3 commercials.



Ordering Information

Component	Cat #	Cat #	Cat #

Contact us

Genei Laboratories Pvt Ltd.

- 📍 No.6, 6th Main
BDA Industrial Suburb
Near SRS Road, Peenya - 560058
- ☎ +91-080 2839 6894
+91-080 2839 1453
- ✉ techsupport@geneilabs.com
info@geneilabs.com
sales@geneilabs.com
www.geneilabs.com